

Meeting Needs. Exceeding Expectations.

630.823.4700 www.iprinc.net

IPR TOOLING INSPECTION & ANALYSIS SYSTEM (TIAS) PRODUCT OVERVIEW

AUTOMATION | UPGRADES | REBUILDS | REPLACEMENT PARTS | FIELD SERVICE | TRAINING | CALIBRATION

Official North American OEM Parts and Service Provider for Manesty Tablet Presses and Coaters

CONTENTS

- Hardware Overview
- Tooling Inspection Fixture
- Tooling Inspection Fixture
- Video Microscope
- Software Overview
- Tooling Inspection
- Tooling Inspection Detail
- Paperless System
- Reporting Functions
- TIAS Components
- TIAS CUSTOMERS
- TIAS Benefits Include
- Price & Availability



630.823.4700 www.iprinc.net

System Overview

Designed For A Tablet Tool Room

- Manage a Complete Tooling Inventory
- Perform Accurate and Repeatable Tablet Tool
 Inspections
- Create A Paperless Tool Room System
- Optimize Tool Selection
- Provide Tool Wear Trending
- Provide A User-Friendly Operator Interface

630.823.4700 www.iprinc.net



Automatic Tooling Inspection Fixture

Streamline Measurements with Improved Accuracy & Efficiency



Working Length

Microscope **Orientation &**

> 630.823.4700 www.iprinc.net



Automatic Tooling Inspection Fixture

Engineered For Precision Measurements



- Linear Slide Bearings
- Hardened Stainless Steel
 Inserts
- Three Point Adjustment For Precise Alignment
- Slide Tables with Digital Readout For Accurate & Easy Cup Depth Setup
- Digimatic Indicators with Air Drives

630.823.4700 www.iprinc.net



Video Microscope

Empower Your Vision Inspection Capabilities



- Computer Screen Viewing
- 200X Magnification
- Proper Tool Orientation
- Punch Face Inspection
- Integrated Fluorescent Light Source

630.823.4700 www.iprinc.net



Software Overview

Easy To Use Software Packed with Functionality For Any Tool Room Operator

- Windows 2000 / XP Professional Compatible
- Microsoft Access Forms
- Borland Interbase Database
- Client Server Architecture
- TIAS Inspection Software (Version 5)
- Stand Alone or Network
- 21 CFR Part 11 Compliant

fprincorporated

630.823.4700 www.iprinc.net

Tooling Inspection

Receive NEW Tooling

New Tool Inspection	- X
2 🛛 🥂 🗐 🗊	View SOPInspect UppersInspect LowersInspect Dies
Details Notes	
Inspection Number:	Upper Punch
Inspection Date: 11/4/2002	Tool Count: 0
Inspected By: KPESCETTO	Tool Size:
Purchase Order:	Bin Location:
Set Number:	Lower Punch
Set Condition:	Tool Count:
Receipt Date:	Tool Size:
Site Location:	Bin Location:
Machines:	Die
Control Number:	Tool Count:
🕹 Upper Drawing: 💽 🐨 🐨	Tool Size:
Revision:	Bin Location:
Lower Drawing:	
Revision:	♣ Supplier:
Record: 1 > > > of 1	

Record Set Number Drawing Number **Tool Quantity Tool Size Tool Location** • View SOPs

• Establish Set Inspection Interval

> 630.823.4700 www.iprinc.net



Tooling Inspection

Record Production Information & Track Usage

		ViewInspectInspectInspectSOPUppersLowersDies
Inspection Number:		Upper Punch
Inspection Date: 11/4/2002		Tool Count: 0
Inspected By: KPESCETTC		or Tool Size:
Batch Number:		
Batch Date:		Lower Punch
🕂 Set Number:		Tool Count: 0
Set Condition:		o Tool Size:
Number of Tablets:	0	
Number of Hours:		Die
Tablet Press:		Tool Count:
Product:		Tool Size: 🗾 💽
Avg. Comp. Force:	(kilonewtons)	I TRALLINATION
Comments:		
Z		

Record
Set Number
Drawing Number
Tool Quantity
Tool Size
Tool Location
View SOPs

• Establish Set Inspection Interval

> 630.823.4700 www.iprinc.net



Tooling Inspection Detail

Individual Tool Inspection Screen

ol Inspection - U	pper Punches		Measurement Method Port	: ACTUAL : 1 (4800)
Standard Me	asurements Tool Specifics Addit	ional Measurements	Comments	X /
Batch Load New Tools	Tool ID:	New Tool	Variance from	Start Manual Measure Start
	Head Outer Diameter:	Measurement:	Spec:	Automatic
	Neck Diameter:			Measure
	Barrel Diameter:			Start Concentricity
	Tip Diameter Major Axis:			Measure
	Tip Diameter Minor Axis:			
	Working Length (calculated):			
	Overall Length:			
	Cup Depth:			
	Concentricity (T.I.R.):			
	Head Hardness:			
	Embossing Height:			
	Head To Barrel Concentricity:			
rd: III	1 I I I I I I I I I I I I I I I I I I I	$\langle - \rangle$	\sim \sim \sim	No. Za

- Measure Upper & Lower Tools
- Measure Dies

lacksquare

- **New Inspections**
- **Usage Inspections**
- Display Measurement Results with Variance From Specification
- Alarm If Out of Specification Event Occurs

630.823.4700 www.iprinc.net



Inventory Look-up

Paperless System

incorporated

🖺 Inventory Look Up					×	
	Tools In Set					
Product Code	Product Name	Set Number	Tablets Per Set	Hours Per Set		
1001-10	10mg Large Tablets	TEST-001	0	0 JIPR, Inc.		
\mathbf{D}	(a)	(a) (a))	(a)		
1001-02	2mg Small Tablets	TEST-003		0 JIPR, Inc.		
\mathbf{D}	(a)	(a) (a)) (a)	(a)		
1001-05	5mg Medium Tablets	TEST-002	0	0 JIPR, Inc.		
\mathbf{D}	(a)	ക്രി ക്ര) (a)	(a)		
					•	
Record: II 1	▶ ▶ ▶ ▶ ★ of 3	•		Þ		
• Products		• Bin Loca	tions	•	Fool Description	
Set Numbe	ers	• Tablets P	roduced	• [Drawing Identific	cation
						630.823.4700
						www.iprinc.net
						www.iprine.net
						of the local data in the local data

Tablet Drawing

Define Critical Tooling Measurements Specific To Your Products

# #	Tablet Drawing Maintenance	e		
	2 🛛		Tablet Drawing cannot be modified as it action record(s).	is being referenced in 1
	Specifications Produc	cts Tolerances		
	Tablet Drawing Num	ber: TEST1-U	Obsolete:	
	Drawing Revis	sion: 0		
	Type:	Upper -	Cup Depth:	0.0550
	Drawing Date:	12/20/2002	Tip Diameter Major Axis:	0.4985
	♣ Supplier:	ECD •	Tip Diameter Minor Axis:	0.3105
	Hob Number:	12345A	Die Inner Diameter Major Axis:	0.5000
	Land:	0.0030	Die Inner Diameter Minor Axis:	0.3125
	Tablet Diameter:	0.5000 × 0.3125	Embossing Letters:	Test1-U
	🕂 Shape:	OVAL	Embossing Height:	0.0090
	Drawing Filena	1me : C:\Documents an	nd Settings\KenP\Desktop\Ariana-Dec2	000.jpg 💿
	Title: 10)mg Large Tablet - Upp	er	
	Comments:			
	Tablet Threshold fo	r Set: 50,000,000	Overload Setting: 12.	00 (kilone w tons)
Re	cord: 🚺 🔳 📘 🕨	▶I ▶* of 6		

- Create Drawing Numbers
- Identify Tip Features
- Record Embossing
- Attach Products
- Define Tablet Production Limit
- Attached Image File

630.823.4700 www.iprinc.net



Reporting Functions

Maintain Detailed Inspection Reports For Quick Reference & Retrieval

Punch and Dies Ins	pection Detail Report				12/30/2002
our Company Name Here					
Inspection Number:	1000009	Upper Tool Size:	B TYPE	Standard B-Type Punch	
Purchase Order Number:	Test1	Lower Tool Size:	BTYPE	Standard B-Type Punch	
Inspection Date:	12/20/2002	Die Tool Size:	1.1875" DIE	1.1875" OD Die	
Set Number:	IPR1	Bin Location:	IPR1 (U) TEST2	(L) TEST3 (D)	
Product Name:	Test	Site Location:	Test1		
Inspection Type:	U	Upper Punch Count:	10		
Inspected By:	TIASMANAGER	Lower Punch Count:	10		
Supplier Name:	Test	Die Count:	0		
Upper Drawing Number:	TEST1-U (0)	Machines:			
Lower Drawing Number:	TEST1-L (0)	Comments:	Repeatability		
Tablet Press:	Repeatability	Condition:	Fair		
Number Of Tablets:	1000000	Lot Date:	12/20/2002		
Average Compression Force:	0	Lot Number:	Repeatability		

Customize
 Report
 Categories

Statistical
 Summaries

 Detailed History of Tooling Inspections

Working Barrel Diameter Cup Depth Head Outer Overall Length Tool ID Diameter Length TEST1-U-01 0.7478 0.0473 0.9985 5.2502 5.2029 TEST1-U-02 0.7479 5.2503 5.2038 0.0465 0.9984 TEST1-U-03 0.7478 0.0471 0.9983 5.2502 5.2031 TEST1-U-04 0.7478 0.0462 0.9982 5.2502 5.2040 TEST1-U-05 0.7478 0.0464 0.9982 5.2502 5.2038 TEST1-U-06 0.7477 <u>0.0465</u> 0.9979 5.2501 5.2036 TEST1-U-07 0.7479 0.0468 0.9986 5.2505 5.2037 0.7477 0.9978 5.2501 5.2038 TEST1-U-08 0.0463 0.7477 5.2502 TEST1-U-09 0.0463 0.9978 5.2039

Pagelof 2



AUTOMATION | UPGRADES | REBUILDS | REPLACEMENT PARTS | FIELD SERVICE | TRAINING | CALIBRATION

630.823.4700 www.iprinc.net

Reporting Functions

Tool Matching Report Optimizes Press Setup To Minimize Weight Variation

Matching Report by Product

12/20/2002

Your Company Name Here

Product: TEST	Uppe	r Set Number	Lower Set Number
		IPR1	IPR1
Upper Tool ID	Upper Working Length	Lower Working Length	Lower Tool ID
TEST1-U-04	5.2040	5.1961	TEST1-L-01
TEST1-U-09	5.2039	5.1980	TEST1-L-05
TEST1-U-08	5.2038	5.1980	TEST1-L-03
TEST1-U-05	5.2038	5.2002	TEST1-L-10
TEST1-U-02	5.2038	5.2004	TEST1-L-08
TEST1-U-07	5.2037	5.2010	TEST1-L-02
TEST1-U-06	5.2036	5.2015	TEST1-L-07
TEST1-U-10	5.2031	5.2016	TEST1-L-09
TEST1-U-03	5.2031	5.2029	TEST1-L-06
TEST1-U-01	5.2029	5.2045	TEST1-L-04
Count 10 Avg	5.2036	5.2004	Count: 10

- Select Tool Sets
- Maximize Weight Control

630.823.4700 www.iprinc.net



Reporting Functions

Variance Report Optimizes Tool Selection For Optimal Press Setup

Tool Variance	e Report					12/30/200
our Company Na	ne Here					
Tool Type: U Set	Number: IPR1		Last Tool-Lev	vel Inspectio	n for Set:	1000009
Measurement: Working	I Length	Average:	5.2039	Variance:	-0.001	to .001
Product: Test						
Tools that Fai						
Tool ID	Set Number		Last Inspection		Measurem	
EST1-U-09	IPR1		100	0009		5.2069

Tools that Passed

Tool ID	Set Number	Last Inspection	Measurement
TEST1-U-01	IPR1	1000009	5.2029
TEST1-U-02	IPR1	1000009	5.2038
TEST1-U-03	IPR1	1000009	5.2031
TEST1-U-04	IPR1	1000009	5.2040
TEST1-U-05	IPR1	1000009	5.2038
TEST1-U-06	IPR1	1000009	5.2036
TEST1-U-07	IPR1	1000009	5.2037
TEST1-U-08	IPR1	1000009	5.2038
TEST1-U-10	IPR1	1000009	5.2031

Total tools that passed: 9



AUTOMATION | UPGRADES | REBUILDS | REPLACEMENT PARTS | FIELD SERVICE | TRAINING | CALIBRATION

User Defined
 Variance
 Tolerances

Out-of-Tolerance Tool Identification

> 630.823.4700 www.iprinc.net

TIAS CUSTOMERS

IPR Customer References

Proctor & Gamble	J & J Puerto Rico
KV	Bristol Myers Squibb PR
McNeil PR	Pfizer Morris Plaines
Apotex	NOW Foods
Abbott Labs - North Chicago	Pfizer PR
Cardinal Health	Purdue Wilson NC
Purdue Canada	Weider Nutrition
Teva	Pfizer Morris Plains
Pfizer Lincoln	Pfizer Brooklyn
Genpharm	Pfizer Canada
Yamanouchi	Pfizer Groton
Wyeth GPR	Schering-Plough
Biovail	Mylan WV
Pfizer Ireland	Pfizer Ann Arbor
Novartis Spain	Abbott Labs - Abbott Park
Alza	Aventis-HMR
Wyeth NY	Barr Labs (Duramed)

Mallinckrodt McNeil Ft. Washington Merck PR Novartis Canada Novartis Lincoln Novartis Suffern Novopharm Perrigo Allegan Perrigo Greenville, SC Searle Brasil Searle Mt. Prospect Searle Skokie Syntex ICN PR Bayer Meyerstown Leiner (Granutec)

630.823.4700 www.iprinc.net



TIAS Benefit Summary

Meeting Needs...Exceeding Expectations!

- Improved Quality Assurance
- Accurate Audit Trail
- Tighter Inventory Control
- Streamlined Inspection Process
 - New Tooling
 - -In-Process Inspections
- Improved Planning, Troubleshooting, and Design
- Reduced Number of Urgent Purchases

630.823.4700 www.iprinc.net

